Identifikačný systém pre enterobaktérie, Identification system for *Enterobacteriaceae*,

SK

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EN**SUMMARY AND EXPLANATION**

ENT 16 is standardized identification system for common species of *Enterobacteriaceae*, it is based on 16-18 miniaturized biochemical tests and software database. List of all microorganisms, for which is kit determined, is placed at the end of leaflet.

PRINCIPLE

Kit consists of 16 wells strip of microtitration plate in classic 96 wells format containing dehydrated substrates, where ENT 16 sp is in the form of strippable microtitration plates and ENT 16 fp is in the form of undivided full plate. Reconstitution of substrates is carried by inoculation of bacterial suspension. During incubation occurs colour change in wells because of metabolic activity of microorganisms. Evaluation of results of tests is visual on base of colour scheme, or by colour formulation described in leaflet. Results of identification can be obtained from evaluation table or by help of evaluation software, which can be found at www.diagnostics.sk/idmicro.

KIT CONTAINS - ACCORDING TO PACKAGING 60 TESTS (sp) / 100 TESTS (fp)

- 10 / 25 microtitration plates of ENT 16
- 60 / 150 result forms
- 10 / 25 incubation packets
- 1 information leaflet

REQUIRED MATERIAL**Reagents:**

- Unbuffered saline 0,85 %, 3,5-5 ml
- Paraffin oil (Ref. 3001)
- IND reagent (Ref. 3002)
- VP and VP reagent (Ref. 2004 and 3004)
- PYR and PYR reagent (Ref. 2003 and 3003)
- Identification software

Materials:

- Pipettes, tampons, loops, burner, tubes and other basic laboratory equipment.

WARNINGS AND SAFETY PRECAUTIONS

- **For in vitro diagnostics use and microbial control**
- **For professional use only**
- Follow the instructions exactly!

- Used strips should be considered as potentially infectious and this must be respected when handling.
- Observe common safety measures according to the regulations of your country.
- Before use, check if the packaging is intact. Do not use damaged kit.

By interpretation of results anamnesis of a patient, source of sample, morphology of colony, microscopic morphology of strain and if necessary, results of all performed tests, mainly results of antibiogram have to be considered.

STORAGE CONDITIONS

Diagnostic kits are delivered in multilayer packets on base of aluminium and organic polymers. Part of each packet is silica gel desiccant. Store kits at temperature from +2 to +25°C. Expiration date is placed on each packaging.

Put unused opened microtitration plate in prepacked Al packet with original silica gel desiccant, close packet carefully and let it by laboratory temperature. Product can be stored two weeks or to date of expiration, in case when it passes earlier.

SAMPLES

Isolate microorganisms, which have to be identified, from suitable unselective cultivation medium (for example blood agar, trypton - soya agar etc.) according to standard microbiological techniques. Make Gram colouring and microscopy of pure culture. Make test of cytochromoxidase – OXI. Identify confirmed isolates on ENT 16 kit.

RECCOMENDED PROCEDURE**Preparation of inoculum**

- Open the tube with saline or use sterile 0,85% solution of sodium chloride in deionized / distilled water.
- Take same well isolated colonies by inoculation loop from pure 18-24 hours old culture.
- Prepare well homogenized suspension of turbidity of 2 McF. Use suspension immediately after preparation.

TIP: Eventually make proof of purity by braying of same loop or tampon, by which was suspension done.

Prepared Petri dish can be used next day to make other additional tests!

Preparation of microtitration plate

- Prepare microtitration plate
- Mark strips with numbers of examined cultures.

Inoculation

- Inoculate 0,1 ml of well homogenized suspension into each well of strip.

Cover tests from URE to LYS (wells from H to D) with three drops of paraffin oil.

Incubation

Put the microtitration plate to packed PE packet, then bend end of packet under plate – this avoid drying of bacterial suspension. Incubate by common atmosphere at temperature $35 \pm 2^\circ\text{C}$ for 18 – 24 hours.

TIP: For optimal incubation conditions, make higher humidity in incubator, for example by adding beaker of clear water or incubate in incubator with managed humidity.

EVALUATION AND INTERPRETATION

Evaluate the tests after incubation with help of interpretation table, colour scheme, or by results of control strain.

Test ONP / IND are bifunctional so after evaluation of primary reaction another result can be obtained from evaluated well of microtitration plate.

Test GLR / ESL is bifunctional also, but there is no need of colorized reagent. It is impossible to get both tests positive.

Make following bifunctional test:

Well H1: ONP / IND - add 2-3 drops of IND reagent and wait 1-2 minutes for colouring of test.

Well A8: GLR / ESL - evaluate positive test GLR for *E. coli*. Mark the negative result of the test GLR also.

- Write down result of bifunctional tests to result form or to evaluate software.

IDENTIFICATION

Result of identification can be obtained by:

- **Identification by identification table:**

Compare results of tests by evaluation table and make identification by results of tests in identification table (page 6).

Identification by identification software:

Enter results of individual tests. When some of the test cannot be evaluated, it can be excluded. Software allows inputting additional tests and so it increases identification effectiveness.

Software microID is available free for our customers on the website of the company or in internet independent form.

QUALITY CONTROL

Quality of produced diagnostic kits is systematically controlled. Chemicals are bought exclusively from certified companies and quality of these chemicals is confirmed by analytical certificate. The functionality of the kits is tested by control collection strains, controlled and tested is also present of bacterial contamination. Kits are exposed to tests of higher temperatures and samples of each batch are saved for right advisement of later reclamations.

For the need of own proof function, use recommended bacterial strains (page 5).

CONSTRAINTS OF METHOD AND MOST OFTEN CAUSES OF WRONG IDENTIFICATION

- Diagnostic kit ENT 16 is determined for identification of bacteria named in this leaflet only.
- Only pure culture of microorganism can be used.
- Tests were not covered by paraffin oil.
- Contamination of wells by inoculum of next strip.
- Used culture is atypical strain.
- Some point of leaflet was not kept.

DETERMINATION CHARACTERISTICS

86 collection strains, strains of clinical background mentioned in database were tested:

88,4 % of tested bacterial strains were identified correctly into rod

84,9 % of tested bacterial strains were identified correctly into species.

WASTE LIQUIDATION

Work with material as with potentially infectious agents. Liquidate leavings according to internal procedures and legislative of your country.

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