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Dátum poslednej revízie: 20.5.2022

SUMMARY AND EXPLANATION

Diagnostic discs VK with content of 5 µg vancomycin and 10 µg colistin are intended for selective cultivation of *Neisseria meningitidis* on base of its resistance to antibiotics contained in disc VK.

PRINCIPLE

Antibiotics diffuse from disc to agar, what creates partly selective zone inhibiting growth of native microbiome.

DIAGNOSTIC KIT CONTAINS

- tube with diagnostic discs - min. 500 pcs (Ref. 6029-500),
- min. 50 pcs (Ref. 6029-50)

WARNINGS AND SECURITY PRECAUTIONS

- **For in vitro diagnostics use.**
- **For professional use only**
- Follow the instructions exactly!
- All samples and inoculated products have to be considered as potentially infectious and obvious safety valid regulations must be respected while handling.
- Do not use product out of expiration.
- Before use, check if the packaging is intact. Do not use damaged kit.
- Leftovers of expired or devaluated product liquidate according to regulations of your country. Packaging materials are recyclable.

STORAGE

Keep discs in dark at temperatures between +2 and +8 °C.

REQUIRED MATERIALS

- Petri dish with blood agar
- loops, burner, tubes and other basic equipment of microbiological laboratory

RECOMMENDED PROCEDURE

1. Spread tested strain on Petri dish.
2. Place sterilely disc VK on inoculated Petri dish, press it gently, so that whole surface of disc is in contact with medium. Incubate for 18-24 h in CO₂ atmosphere at temperature 35 ± 2 °C.

EVALUATION AND INTERPRETATION


In inhibition zone can grow yeasts (differentiation by microscopy), *Neisseria lactamica* (ONP test positive) and *Neisseria meningitidis* (confirm by NEISS 8)

Functionality of diagnostic discs can be verified by these control strains:

Bacterial strain	CCM / ATCC	Inhibition zone
<i>Neisseria meningitidis</i>	4634 / 13077	no zone
<i>Streptococcus pneumoniae</i>	4501 / 49619	> 20 mm

CCM: Czech Collection of Microorganisms, Masaryk University Brno, Kamenice 5, 625 00 Brno, Czech Republic, tel. +420549491430, e-mail: ccm@sci.muni.cz

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