

## PP-20/4, platform

### DESCRIPTION

Flat platform with non-slip rubber mat can accommodate various low profile containers.

### SPECIFICATIONS

|                            |                    |
|----------------------------|--------------------|
| Working area               | 355 × 455 mm × 4x  |
| Height between platforms   | 140 mm             |
| Overall dimensions (W×D×H) | 380 × 480 × 440 mm |



### CAT. NUMBER

BS-010126-EK

Platform PP-20-4



PSU-20i  
Without platform  
Multi-functional Orbital Shaker

Shaker PSU-20i provides three motion types: orbital, reciprocal and vibrating, which can be performed separately, pairwise and sequentially in repeated cycles. Shaker is designed for applications both in small specialized ...

[read more](#)