

PRSC-10, platform

DESCRIPTION

10 tubes 25-30 mm diameter (50 ml tubes).

SPECIFICATIONS

| | |
|---------------|----------|
| Platform type | optional |
| Tube diameter | 20-30 mm |
| Tube volume | 50 ml |
| Capacity | 10 |



CAT. NUMBER

BS-010117-JK

Platform PRSC-10



Multi Bio RS-24
Including PRS-26 platform
Programmable rotator

Programmable Rotator **Multi Bio RS-24** performs several motion types in one module. Microprocessor control allows performing not only Vertical rotation (1) of the platform, but also Reciprocal rotation (2) as ...

[read more](#)