

# PRSC-18, platform

## DESCRIPTION

Can accommodate 18 tubes, 16 mm diameter (15 ml).

## SPECIFICATIONS

Platform type	optional
Tube diameter	16 mm
Tube volume	15 ml
Capacity	18



## CAT. NUMBER

BS-010117-EK

Platform PRSC-18



**Bio RS-24**  
Including PRS-22 platform  
Mini-rotator

Mini-rotator **Bio RS-24** provides vertical rotation of the platform. The rotator is an ideal instrument for preventing blood coagulation in tubes and for fulfilment of procedures of biological components extraction. ...

[read more](#)