

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321 13088 Berlin, Germany

T: +49 30 700 144-0 F: +49 30 700 144-30 info@sifin.de www.sifin.de

Product Anti-Coli O 126:K 71

Catalogue no. TS2212

Package size 1 ml

Lot no. 9141122

Expiry date 2027-10

Intended use Anti-Coli O 126:K 71 records all *E. coli* strains that contain the antigen

O 126:K 71 in slide agglutination. In the Widal reaction the serum reaches a titre of 1:1600 (tube and MTP) with the homologous *E. coli* strains; deviations

of ±1 titre step are acceptable.

Specificity Since unspecific co-agglutinations with other E. coli strains often occur in the

slide agglutination test, a positive result in the slide test must always be confirmed by examination of the boiled *E. coli* isolate in the Widal test.

Preservatives Sodium azide 0.9 mg/ ml

Storage 2 °C...8 °C