

# Liofilchem®

## Certificate of Analysis

Page 1 of 1

<b>Product</b>	<b>Batch</b>	<b>Expiration date</b>
Sulfonamide S3 300 µg	070517106	2020.07.04
<b>Ref.</b>	<b>Initial</b>	
9126 – 9126/1	S3	

### Antimicrobial Susceptibility Testing

Tested according to current CLSI methodology

<b>Control strains</b>	<b>Medium</b>	<b>Inoculum</b>	<b>Incubation</b>	<b>Expected Results</b> Zone range (mm)	<b>Results</b> Zone (mm)
<i>Escherichia coli</i> ATCC® 25922	Mueller Hinton II Agar	0.5 McFarland in saline	35 ± 2°C, ambient 16-18 h	15-23	23
<i>Staphylococcus aureus</i> ATCC® 25923	Mueller Hinton II Agar	0.5 McFarland in saline	35 ± 2°C, ambient 16-18 h	24-34	29

### Batch Release

Approved

**Date**

10.07.2017

**Signature**Quality Control  
(D. Vitagliano)

The results reported were obtained at the time of release.