

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
 ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321  
 13088 Berlin, Germany

T: +49 30 927 030-0  
 F: +49 30 927 030-30  
 info@sifin.de  
 www.sifin.de

<b>Product</b>	<b>Anti-Salmonella O:46</b>
Catalogue no.	TR 1315
Package size	1 ml
<b>Lot no.</b>	<b>191 10 14</b>
<b>Expiry date</b>	<b>2018-10</b>
Intended use	Anti-Salmonella O:46 records all <i>Salmonella</i> strains of group D <sub>2</sub> (O:9,46) in the slide agglutination test. Strains of group D <sub>3</sub> (O:9,12,46,27) do not agglutinate or agglutinate only weakly.
Specificity	No agglutination with heterologous <i>Salmonella</i> O antigens and <i>Salmonella</i> H antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,  
 place, date

10.11.14

signature Head of Quality Control

