

Certificate of Quality Control

sifin diagnostics gmbh

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

Berliner Allee 317-321
13088 Berlin, Germany

T: +49 30 927 030-0
F: +49 30 927 030-30
info@sifin.de
www.sifin.de

Product	Anti-Salmonella O:46
Catalogue no.	TR 1315
Package size	1 ml
Lot no.	0270118
Expiry date	2022-01
Intended use	Anti-Salmonella O:46 records all <i>Salmonella</i> strains of group D ₂ (O:9,46) in the slide agglutination test. Strains of group D ₃ (O:9,12,46,27) do not agglutinate or agglutinate only weakly.
Specificity	No agglutination with heterologous <i>Salmonella</i> O antigens and <i>Salmonella</i> H antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,
place, date

5.02.2018


signature Head of Quality Control