

Certificate of Quality Control

Manufacturer **SIFIN** Institut für Immunpräparate und Nährmedien GmbH Berlin,
Berliner Allee 317-321, 13088 Berlin, Germany

Product **Enteroclon Anti-Salmonella O 35**

Catalogue No. TR 1341 (1 ml)
(Package size)

Lot No. **68 06 13**

Expiry date **2018-06**

Intended use Enteroclon Anti-Salmonella O 35 records all Salmonella strains, as long as they contain antigen O 35.

Specificity No other agglutinations with heterologous Salmonella O antigens and Salmonella H antigens are known.

Preservatives Sodium azide 0.9 mg / ml

Storage 2...8 °C

Berlin,

01.07.13


Head of Quality Control