

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
 ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
 13088 Berlin, Germany

T: +49 30 927 030-0
 F: +49 30 927 030-30
 info@sifin.de
 www.sifin.de

Product	Anti-Shigella sonnei F-form (phase II)
Catalogue no.	TR 2013
Package size	1 ml
Lot no.	7370817
Expiry date	2022-08
Intended use	Anti-Shigella sonnei F-form (phase II) records the strains of <i>Shigella sonnei</i> F- form (phase II) in slide agglutination.
Specificity	No agglutination with antigens of <i>Shigella sonnei</i> S-form (phase I), <i>Shigella flexneri</i> , <i>Shigella dysenteriae</i> (type 1 to 10) and <i>Shigella boydii</i> (type 1 to 15).
Preservatives	Sodium azide 0.9 mg / ml
Storage	2 °C...8 °C

Berlin,

25.09.2017

place, date



signature Head of Quality Control